


**Thermo Scientific**  
B·R·A·H·M·S MR-proADM



## MR-proADM

# Optimize risk assessment, decision making and treatment processes

High accuracy in predicting short- and long-term outcomes

- Established cut-offs
- Improved patient classification

Better alignment of clinical resources

- Excellent prognosis profile in sepsis
- Used in a wide range of clinical applications

**Thermo**  
SCIENTIFIC



# Optimize risk assessment, decision making and treatment processes

MR-proADM is a blood biomarker that provides accurate short-, mid- and long-term prognostic information and aids in the triage and multi-dimensional risk assessment of patients in the Emergency Department (ED) and Intensive Care (ICU) settings.


## Adrenomedullin is the key

Adrenomedullin, a 52 amino acid peptide, is a member of the calcitonin peptide family<sup>1</sup> and is widely expressed in many tissues and organs. It has been shown to have a variety of physiological functions, including immune-modulating, direct bactericidal, diuretic and potent vasodilatory activity<sup>1-5</sup>, and in healthy conditions, circulates at low picomolar concentrations<sup>6</sup>.

In many diseased states such as hypertension<sup>7</sup>, renal failure<sup>7</sup>, lower respiratory diseases<sup>8</sup> and septic shock<sup>9</sup>,

plasma levels are significantly up-regulated in proportion to disease severity<sup>7, 10-13</sup>. This allows clinicians to determine the patients **most at risk of developing complications** on admission or during ED and ICU stay, in order to rapidly triage and **administer the most effective treatment**, in the shortest possible time.

**The clinical and financial benefits of an early diagnosis and risk assessment to both the patient and healthcare provider, therefore, cannot be overstated.**





# Our solution: MR-proADM

## Fast, stable and independent



High correlation with increasing severity of disease<sup>17</sup>, and accurate short- to long-term prognostic prediction capability in the ED and ICU

MR-proADM is a novel biomarker which provides more precise patient risk management and greater confidence in treatment site assignment.

Used in a variety of indications, it can enhance routine clinical investigation and treatment, and provide a viable alternative to many current risk assessment scores.

### Use in the Emergency Department (ED)

Rapid triaging and risk assessment of patients on admission and throughout the ED can:

- Decrease time to treatment
- Increase out-patient numbers
- Reduce length of stay

### Use in the Intensive Care Unit (ICU)

Immediate risk assessment of mortality and adverse effects in the ICU can:

- Maximize patient safety
- Guide the most appropriate treatment
- Provide an early warning of additional complications

## Optimization of risk assessment and patient management using MR-proADM



- Levels are **not** influenced by food or water intake<sup>18</sup>



- No significant **gender** related differences<sup>18</sup>



- **Stability** of up to 72 hours in EDTA plasma at room temperature<sup>18</sup> and over four freeze/thaw cycles<sup>18</sup>



- Well **documented** for clinical use



- **Rapidly** available to aid timely clinical decision-making using the KRYPTOR™ platform\*

\* assay incubation time 29 mins and small sample volume (26 µL in EDTA plasma)



# Mechanism of action and functions of Adrenomedullin



<div>DOWN</div> <div>Adrenomedullin is down-regulated by</div> <div><div>IFN-<math>\gamma</math></div><div>TGF-<math>\beta</math></div></div>	<div>DOWN</div> <div>Adrenomedullin down-regulates</div> <div><div>TNF-<math>\alpha</math></div><div>IL-6</div><div>IL-1<math>\beta</math></div><div>IL-12</div><div>IFN-<math>\gamma</math></div></div>	<div>UP</div> <div>Adrenomedullin is up-regulated by</div> <div><div>LPS</div><div>TNF-<math>\alpha</math></div><div>IL-6</div><div>Hypoxia</div></div>	<div>UP</div> <div>Adrenomedullin up-regulates</div> <div><div>Shear stress</div><div>IL-1<math>\beta</math></div><div>Nitric oxide</div><div>IL-10</div></div>
---	--	---	---

**Table 1:**  
Molecules which regulate and are regulated by ADM.


Adrenomedullin, present mainly in endothelial and vascular smooth muscle cells<sup>19</sup>, can act as both a hormone and cytokine (often termed a “hormokine”<sup>20</sup>) in an autocrine and paracrine manner<sup>21</sup>. Its potent vasodilatory and hypotensive response is elicited through an initial increase in cyclic adenosine monophosphate levels, and a subsequent production of nitric oxide<sup>1, 22, 23</sup>.

The importance of adrenomedullin in homeostasis is illustrated by its central role in the up- and down-regulation of cytokines and other mediators (table 1), as well as its own stimulatory and inhibitory<sup>24, 25</sup> effect on cytokine production. Indeed, Interleukin (IL)-1 $\beta$  and tumour necrosis factor (TNF) are two of the most potent stimulators for adrenomedullin production<sup>26</sup> and adrenomedullin itself is up-regulated by hypoxia, bacterial products and shear stress<sup>20</sup>.

The ubiquitous and important functional role of adrenomedullin results in its clinical use in many diverse indications (figure 2). For example, its precursor molecule, MR-proADM, has been shown to be a powerful risk assessment marker in sepsis and lower respiratory tract infections<sup>8, 17, 27, 28</sup>, with the ability to predict 30 day mortality regardless of the underlying diagnosis<sup>29</sup>, and in the non-specific complaints of elderly patients<sup>30</sup>. Its plasma concentrations have been shown to be elevated in myocardial infarction, and to correlate with the severity of acute and chronic heart failure<sup>31-34</sup>, as well as being elevated in various types of glomerulonephritis and progressively increased in patients with chronic renal failure<sup>35, 36</sup>.



# MR-proADM: Used in a wide range of clinical applications



**Figure 2:**  
Tissue specific functions of Adrenomedullin and different indications applicable to MR-proADM.



# Cut-off values: For more accurate risk assessment and prognosis



Established cut-off values are important to help guide a clinician into making a correct assessment of risk. Healthy individuals have detectable levels of MR-proADM of approximately 0.4 nmol/L, but this value can increase

significantly, depending on individual disease conditions. Cut-off values for conditions such as severe sepsis, cardiogenic shock and COPD can be found in figure 3, along with corresponding references.




Figure 3: Published cut-off levels for MR-proADM in different indications.



# Faster and more accurate patient risk stratification



The combination of established MR-proADM cut-off values with current risk assessment scores offers an easy to use algorithm for the more accurate assessment of a patient's risk of clinical complications, and a greater confidence in determining the most appropriate site of care (figure 4).

Using such an approach, **significantly more patients can be classified as having a low risk of future complications**, leading to an increase in those that can be discharged safely and treated as out-patients<sup>37</sup>. This consequently reduces the financial and resource consumption pressure on the healthcare provider, whilst minimizing any unnecessary clinical complications for the patient through an unwarranted hospital stay.

Furthermore, the addition of MR-proADM not only **more accurately highlights patients who should be classified as an intermediate risk**, through either a decrease or increase in risk severity; but it can also **more accurately highlight those who are at greatest risk of developing complications**, either on admission or throughout their hospital stay. Indeed, studies have shown

that using MR-proADM, significant numbers of patients, formerly classified as high risk, can be treated at a lower intensity or even as out-patients<sup>30, 37</sup>, thus freeing up precious resources in either the ICU or other high risk settings.

Accordingly, more accurate patient classification due to the addition of MR-proADM can lead to<sup>30, 38</sup>:

- Increased out-patient numbers
- Enhanced discharge management
- Avoiding under-treatment of patients at risk
- Decreased length of stay
- Reduced adverse complications
- Reduced costs

		1. Initial Clinical Risk Assessment by ED Physician (ITT)		
		Low Risk	Intermediate Risk	High Risk
2. MR-proADM	Low ( $< 0.75$ nmol/L)	Low	Low	Intermediate
	Intermediate ( $\geq 0.75$ and $\leq 1.5$ nmol/L)	Low	Intermediate	High
	High ( $> 1.5$ nmol/L)	Intermediate	High	High

**Figure 4:** Recommended algorithm for MR-proADM in combination with established clinical risk assessment scores in patients with, for example, non-specific complaints<sup>39</sup>.



## References

- Kitamura K, Kangawa K, Kawamoto M, et al. Adrenomedullin: a novel hypotensive peptide isolated from human pheochromocytoma. *Biochem Biophys Res Commun.* Apr 30 1993;192(2):553-560.
- Hinson JP, Kapas S, Smith DM. Adrenomedullin, a multifunctional regulatory peptide. *Endocr Rev.* Apr 2000;21(2):138-167.
- Linscheid P, Seboek D, Zulewski H, Keller U, Muller B. Autocrine/paracrine role of inflammation-mediated calcitonin gene-related peptide and adrenomedullin expression in human adipose tissue. *Endocrinology.* Jun 2005;146(6):2699-2708.
- Nishikimi T. Adrenomedullin in the kidney-renal physiological and pathophysiological roles. *Curr Med Chem.* 2007;14(15):1689-1699.
- Pio R, Martinez A, Unsworth EJ, et al. Complement factor H is a serum-binding protein for adrenomedullin, and the resulting complex modulates the bioactivities of both partners. *J Biol Chem.* Apr 13 2001;276(15):12292-12300.
- Lim SC, Morgenthaler NG, Subramaniam T, Wu YS, Goh SK, Sum CF. The relationship between adrenomedullin, metabolic factors, and vascular function in individuals with type 2 diabetes. *Diabetes Care.* Jun 2007;30(6):1513-1519.
- Ishimitsu T, Nishikimi T, Saito Y, et al. Plasma levels of adrenomedullin, a newly identified hypotensive peptide, in patients with hypertension and renal failure. *J Clin Invest.* Nov 1994;94(5):2158-2161.
- Schuetz P, Wolbers M, Christ-Crain M, et al. Prohormones for prediction of adverse medical outcome in community-acquired pneumonia and lower respiratory tract infections. *Crit Care.* 2010;14(3):R106.
- Nishio K, Akai Y, Murao Y, et al. Increased plasma concentrations of adrenomedullin correlate with relaxation of vascular tone in patients with septic shock. *Crit Care Med.* Jun 1997;25(6):953-957.
- Nishikimi T, Saito Y, Kitamura K, et al. Increased plasma levels of adrenomedullin in patients with heart failure. *J Am Coll Cardiol.* Nov 15 1995;26(6):1424-1431.
- Washimine H, Yamamoto Y, Kitamura K, et al. Plasma concentration of human adrenomedullin in patients on hemodialysis. *Clin Nephrol.* Dec 1995;44(6):389-393.
- Kobayashi K, Kitamura K, Hirayama N, et al. Increased plasma adrenomedullin in acute myocardial infarction. *Am Heart J.* Apr 1996;131(4):676-680.
- Kakishita M, Nishikimi T, Okano Y, et al. Increased plasma levels of adrenomedullin in patients with pulmonary hypertension. *Clin Sci (Lond).* Jan 1999;96(1):33-39.
- Meeran K, O'Shea D, Upton PD, et al. Circulating adrenomedullin does not regulate systemic blood pressure but increases plasma prolactin after intravenous infusion in humans: a pharmacokinetic study. *J Clin Endocrinol Metab.* Jan 1997;82(1):95-100.
- Lewis LK, Smith MW, Brennan SO, Yandle TG, Richards AM, Nicholls MG. Degradation of human adrenomedullin(1-52) by plasma membrane enzymes and identification of metabolites. *Peptides.* 1997;18(5):733-739.
- Struck J, Tao C, Morgenthaler NG, Bergmann A. Identification of an Adrenomedullin precursor fragment in plasma of sepsis patients. *Peptides.* Aug 2004;25(8):1369-1372.
- Christ-Crain M, Morgenthaler NG, Struck J, Harbarth S, Bergmann A, Muller B. Mid-regional pro-adrenomedullin as a prognostic marker in sepsis: an observational study. *Crit Care.* 2005;9(6):R816-824.
- Morgenthaler NG, Struck J, Alonso C, Bergmann A. Measurement of midregional proadrenomedullin in plasma with an immunoluminometric assay. *Clin Chem.* Oct 2005;51(10):1823-1829.
- Kitamura K, Kangawa K, Eto T. Adrenomedullin and PAMP: discovery, structures, and cardiovascular functions. *Microsc Res Tech.* Apr 1 2002;57(1):3-13.
- Schuetz P, Marlowe RJ, Mueller B. The prognostic blood biomarker proadrenomedullin for outcome prediction in patients with chronic obstructive pulmonary disease (COPD): a qualitative clinical review. *Clin Chem Lab Med.* Sep 25 2014.
- Elasser TH, Kahl S. Adrenomedullin has multiple roles in disease stress: development and remission of the inflammatory response. *Microsc Res Tech.* Apr 15 2002;57(2):120-129.
- Shimekake Y, Nagata K, Ohta S, et al. Adrenomedullin stimulates two signal transduction pathways, cAMP accumulation and Ca<sup>2+</sup> mobilization, in bovine aortic endothelial cells. *J Biol Chem.* Mar 3 1995;270(9):4412-4417.
- Eto T. A review of the biological properties and clinical implications of adrenomedullin and proadrenomedullin N-terminal 20 peptide (PAMP), hypotensive and vasodilating peptides. *Peptides.* Nov 2001;22(11):1693-1711.
- Kubo A, Minamino N, Isumi Y, et al. Production of adrenomedullin in macrophage cell line and peritoneal macrophage. *J Biol Chem.* Jul 3 1998;273(27):16730-16738.
- Wong LY, Cheung BM, Li YY, Tang F. Adrenomedullin is both proinflammatory and antiinflammatory: its effects on gene expression and secretion of cytokines and macrophage migration inhibitory factor in NR8383 macrophage cell line. *Endocrinology.* Mar 2005;146(3):1321-1327.
- Sugo S, Minamino N, Shoji H, et al. Interleukin-1, tumor necrosis factor and lipopolysaccharide additively stimulate production of adrenomedullin in vascular smooth muscle cells. *Biochem Biophys Res Commun.* Feb 6 1995;207(1):25-32.
- Abraham E, Laterre PF, Garg R, et al. Drotrecogin alfa (activated) for adults with severe sepsis and a low risk of death. *N Engl J Med.* Sep 29 2005;353(13):1332-1341.
- Huang DT, Angus DC, Keillum JA, et al. Midregional proadrenomedullin as a prognostic tool in community-acquired pneumonia. *Chest.* Sep 2009;136(3):823-831.
- Potocki M, Br idhardt T, Reichlin T, et al. Midregional pro-adrenomedullin in addition to b-type natriuretic peptides in the risk stratification of patients with acute dyspnea: an observational study. *Crit Care.* 2009;13(4):R122.
- Nickel CH, Messmer AS, Geigy N, et al. Stress markers predict mortality in patients with nonspecific complaints presenting to the emergency department and may be a useful risk stratification tool to support disposition planning. *Acad Emerg Med.* Jul 2013;20(7):670-679.
- Klip IT, Voors AA, Anker SD, et al. Prognostic value of mid-regional pro-adrenomedullin in patients with heart failure after an acute myocardial infarction. *Heart.* Jun 2011;97(11):892-898.
- Travaglino F, Russo V, De Berardinis B, et al. Thirty and ninety days mortality predictive value of admission and in-hospital procalcitonin and mid-regional pro-adrenomedullin testing in patients with dyspnea. Results from the VErifying DYspnea trial. *Am J Emerg Med.* Apr 2014;32(4):334-341.
- von Haehling S, Filippatos GS, Pappasotiropoulos J, et al. Mid-regional pro-adrenomedullin as a novel predictor of mortality in patients with chronic heart failure. *Eur J Heart Fail.* May 2010;12(5):484-491.
- Haaf P, Twerenbold R, Reichlin T, et al. Mid-regional pro-adrenomedullin in the early evaluation of acute chest pain patients. *Int J Cardiol.* Sep 30 2013;168(2):1048-1055.
- Ishihara T, Yokota N, Hisanaga S, et al. Increased plasma levels of mature form of adrenomedullin in patients with chronic renal failure. *Clin Nephrol.* Aug 1999;52(2):119-123.
- Kinoshita H, Fujimoto S, Kitamura K, et al. Increased plasma levels of mature adrenomedullin in chronic glomerulonephritis. *Nephron.* Nov 2000;86(3):333-338.
- Albrich WC, Dusemund F, Ruegger K, et al. Enhancement of CURB65 score with proadrenomedullin (CURB65-A) for outcome prediction in lower respiratory tract infections: derivation of a clinical algorithm. *BMC Infect Dis.* 2011;11:12.
- Albrich WC, Ruegger K, Dusemund F, et al. Optimised patient transfer using an innovative multidisciplinary assessment in Kanton Aargau (OPTIMA I): an observational survey in lower respiratory tract infections. *Swiss Med Wkly.* 2011;141:w13237.

To request more information, please email: [info.proadm@thermofisher.com](mailto:info.proadm@thermofisher.com) or contact your local account manager

**thermoscientific.com**

© 2015 Thermo Fisher Scientific Inc. All rights reserved.

KRYPTOR is a registered trademark of CIS bio international, licensed for use by B-R-A-H-M-S, a part of Thermo Fisher Scientific. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. All data regarding specifications, terms and pricing correspond to the existing knowledge at the time of the printing. We are not responsible for any errors, misprints or changes. Reprint, also in parts, solely with prior written consent of B-R-A-H-M-S GmbH.

Thermo Fisher Scientific products are distributed worldwide; not all are intended uses and applications mentioned in this printing are registered in every country.

The manufacture and/or use of this product is covered by one or more of the following patents: EP1488209, EP2347266, MX255044, JP4602321, CN200480006345, EP2097748

**Clinical Diagnostics**

Thermo Fisher Scientific  
B-R-A-H-M-S GmbH  
Neuendorfstr. 25  
16761 Hennigsdorf  
Deutschland

[www.thermoscientific.com/brahms](http://www.thermoscientific.com/brahms)  
[www.thermoscientific.com/proadrenomedullin](http://www.thermoscientific.com/proadrenomedullin)

+49 (0)3302 883 0  
+49 (0)3302 883 100 Fax  
[info.brahms@thermofisher.com](mailto:info.brahms@thermofisher.com)

**Thermo**  
SCIENTIFIC

A Thermo Fisher Scientific Brand