



# concile® Ω100

fast, simple, portable.

- rapid and instant point of care testing (POCT)
- easy to use, calibration-free
- mobile use in clinics, hospitals and ambulances

# concile® Ω100

The reader for easy and accurate quantitative point of care testing (POCT).

## 1. Enter patient data



## 2. Apply sample



## 3. Insert the test device



## 4. Automatic measurement after 15 minutes



## 5. Read value, print, save



## Benefits at a glance

- Automatic identification of the test parameters via barcode
- Calibration-free: Calibration data on SD card supplied
- In accordance to RiliBÄK: „The 2008 Directive of the German Medical Association on the Quality Assurance of Tests in Laboratory Medicine“
- Controls available for each parameter
- Computer connectivity via USB cable and software *Patient Data Viewer*
- Network connectivity via HL7 and GDT protocol
- Internal test of devices functioning
- Shock resistant measurement process
- Mobile: Can be operated with environmentally friendly rechargeable batteries
- Comes with printer, software and carrying case

concile Ω100	Specifications
<b>Weight</b>	1,76 lb (800 g)
<b>Dimensions</b>	8.27 x 5.12 x 2.76 inches (210 x 130 x 70 mm)
<b>Display</b>	4.3 inches x TFT color display
<b>Battery</b>	Li-Ion battery, operating time 8 hours
<b>Accessories</b>	Thermal printer accessories Barcode scanner integrated
<b>Memory</b>	1 MB of memory sufficient for 3500 data sets
<b>Connections</b>	RS323 serial interface (printer) USB (computer connection) SD card slot (calibration) LAN port
<b>Menu languages</b>	German, English, Spanish, Italian, French, Portuguese, Swedish, Romanian, Serbo-Croatian

concile GmbH  
Kronenmattenstr. 6  
D-79100 Freiburg i. Brsg.  
T +49 (0)761.15 14 74-0  
F +49 (0)761.15 14 74-19  
info@concile.de  
www.concile.de

concile®  
innovation meets diagnostics