

Thermo Scientific B·R·A·H·M·S Tumor Markers

Tumor marker combinations for therapy control and follow-up*


Marker Tumor	CEA	AFP	CA 19-9	CA 125 II	CA 15-3	Chromo- granin A	NSE	CYFRA 21-1	SCC	hCG+β	PSA	Calci- tonin	hTG
Colon	●		●			●							
Pancreas	●		●			●							
Stomach	●		●										
Esophagus	●		●						●				
Liver (HCC)	○	●											
Biliary ducts			●										
Neuro- endocrine tumors (NET)						●	●						
Breast	●				●								
Ovary				●									
Cervix	●								●				
Chorion										●			
Lung SCLC						○	●	●					
NSCLC	●							●	●				
Germ cell		●								●			
Prostate	○										●		
Bladder	○							●					
Thyroid C-Cell	●					○						●	●
ENT	●							●	●				

● 1st choice ● 2nd choice ○ 3rd choice

*by courtesy of Dr. P. Stieber

Thermo Scientific Tumor Markers

A broad range of markers available



Brain
NSE, CEA
Pituitary Gland
Prolactin
Lung
SCLC: NSE, CYFRA 21-1, Chromogranin A
NSCLC: SCC, CYFRA 21-1, CEA
Liver
AFP, CEA, CA 19-9, Chromogranin A
Gallbladder
CA 19-9, CEA
Adrenal Gland
Chromogranin A
Kidney
CEA, NSE
Colon
CEA, CA 19-9, Chromogranin A
Bladder
CYFRA 21-1, CEA, NSE
Prostate
Total PSA, Free PSA
Testicle
AFP, hCG+ β
Bone Metastases
Osteocalcin

ENT (Ear, Nose, and Throat)
SCC, CYFRA 21-1, CEA
Thyroid
MTC: Calcitonin, Chromogranin A, CEA
DTC: Thyroglobulin, CEA
Esophagus
SCC, CYFRA 21-1, CEA
Breast
CA 15-3, CEA
Stomach
CEA, CA 19-9, Chromogranin A
Pancreas
CA 19-9, CEA, Chromogranin A
Neuroendocrine Tumors
Chromogranin A, NSE
Ovary
CA 125 II, CEA, AFP, hCG+ β
Uterus
SCC, hCG+ β , CYFRA 21-1
Cervix
SCC, CYFRA 21-1, CEA
Vagina, vulva
SCC

thermoscientific.com/brahms

© 2013 Thermo Fisher Scientific Inc. All rights reserved. B-R-A-H-M-S™ and all other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. All data regarding specifications, terms and pricing correspond to the existing knowledge at the time of the printing. We are not responsible for any errors, misprints or changes. Reprint, also in parts, solely with prior written consent of B-R-A-H-M-S GmbH.

Thermo Fisher Scientific products are distributed worldwide; not all intended uses and applications mentioned in this printing are registered in every country.

Clinical Diagnostics

Thermo Fisher Scientific
B-R-A-H-M-S GmbH
Neuendorferstr. 25
16761 Hennigsdorf
Germany

106703.3

+49 (0)3302 883 0
+49 (0)3302 883 100 fax
info.brahms@thermofisher.com

www.thermoscientific.com/brahms
www.thermoscientific.com/copeptin
www.thermoscientific.com/proadrenomedullin
www.thermoscientific.com/procalcitonin
www.thermoscientific.com/kryptor

Thermo
S C I E N T I F I C

Part of Thermo Fisher Scientific